

PRODUCT INFORMATION

Common Name	CIM 331,CIM331
Conjugate	PE-conjugated
Synonyms	CRL; GPL; CRL3; GLMR; GLM-R; PLCA2; hGLM-R; IL-31RA; PRO21384; zcytoR17
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG2(C135S,R137K,E141G,S142G,C223S,H268Q,R355Q,S378A,Q419E) - Kappa
Reactivity	Human
Target	IL31RA
Uniprot ID	Q8NI17
Description	PE-conjugated Anti-IL31RA(nemolizumab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

