

PRODUCT INFORMATION

Common Name	TBC
Conjugate	PE-conjugated
Synonyms	hK2;hGK-1;KLK2A2
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	1 μ l/test
Formulation & Reconstitution	Liquid,PBS with,0.18%BSA□0.1%Proclin300
Host Species	N/A
IgG type	IgG1
Reactivity	Human
Target	KLK2
Uniprot ID	P20151
Description	PE-conjugated Anti-KLK2(alpitatug biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

