

**PRODUCT INFORMATION**

<b>Common Name</b>	ABC-008,ABC008
<b>Conjugate</b>	PE-conjugated
<b>Synonyms</b>	2F1;CLEC15A;MAFA;MAFA-2F1;MAFA-L;MAFA-LIKE
<b>Applications</b>	Flow Cyt
<b>Endotoxin</b>	Less than 1.0 EU/ $\mu$ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Humanized
<b>IgG type</b>	Human IgG1(E356D,M358L) - kappa
<b>Reactivity</b>	Human
<b>Target</b>	KLRG1
<b>Uniprot ID</b>	Q96E93
<b>Description</b>	PE-conjugated Anti-KLRG1(ulviprubart biosimilar) mAb
<b>Delivery</b>	Under Development
<b>Storage&amp;Shipping</b>	Store at 2°C-8°C for 6 months
<b>Sterility</b>	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 $\mu$ m) prior to use.
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only

