

PRODUCT INFORMATION

Common Name	REGN4461
Conjugate	PE-conjugated
Synonyms	LEP-R; CD295; DB; OBR
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens Humanized
IgG type	Human IgG4 - Kappa
Reactivity	Human
Target	LEPR
Uniprot ID	P48357
Description	PE-conjugated Anti-LEPR(mibavademab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

