Cat. No. DMC100277P



PRODUCT INFORMATION

Clone ID DMC277 LGALS1 **Target Synonyms** GAL1; GBP **Host Species** Rabbit

PE-conjugated Anti-LGALS1 antibody(DMC277); Description

IgG1 Chimeric mAb

Delivery Under Development

Uniprot ID P09382

Rabbit/Human Fc chimeric IgG1 IgG type

Clonality Monoclonal Reactivity Human **Applications** Flow Cyt

Recommended Flow Cyt 1:100 **Dilutions**

Background

Purified from cell culture supernatant by affinity **Purification**

chromatography

Formulation & Liquid PBS with 0.05% Proclin300, 1% BSA Reconstitution

Storage & Shipping Store at 2°C-8°C for 6 months

> The galectins are a family of beta-galactosidebinding proteins implicated in modulating cell-cell and cell-matrix interactions. This gene product

may act as an autocrine negative growth factor that regulates cell proliferation.

Research use only **Usage** Conjugate PE-conjugated

> All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or

DIMA Disclaimer reverse engineering attempt is prohibited. We are

actively scrutinizing all patent application to ensure no IP infringement.





