

PRODUCT INFORMATION

Common Name	enfortumab
Conjugate	PE-conjugated
Synonyms	EDSS1;LNIR;nectin-4;PRR4;PVRL4
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 0.5µl/test
Formulation & Reconstitution	Liquid,PBS with,0.18%BSA□0.1%Proclin300
Host Species	Homo sapiens
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	Nectin-4
Uniprot ID	Q96NY8
Description	PE-conjugated Anti-Nectin-4(enfortumab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 0.5µl/test PE-conjugated Anti-Nectin-4(enfortumab biosimilar) on T-47D cell line (Blue histogram) or CHO-S cell line (Red histogram).

