

**PRODUCT INFORMATION**

<b>Common Name</b>	PF-04950615, RN-316, RN316
<b>Conjugate</b>	PE-conjugated
<b>Synonyms</b>	FH3;FHCL3;HCHOLA3;LDLCQ1;NARC-1;NARC1;PC9
<b>Applications</b>	Flow Cyt
<b>Endotoxin</b>	Less than 1.0 EU/ $\mu$ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Humanized
<b>IgG type</b>	Human IgG2 - Kappa
<b>Reactivity</b>	Human
<b>Target</b>	PCSK9
<b>Uniprot ID</b>	Q8NBP7
<b>Description</b>	PE-conjugated Anti-PCSK9(bococizumab biosimilar) mAb
<b>Delivery</b>	Under Development
<b>Storage&amp;Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only

