

PRODUCT INFORMATION

Common Name	MK-3475
Conjugate	PE-conjugated
Synonyms	PDCD1;PD1;CD279;SLEB2
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG4
Reactivity	Human
Target	PD-1
Uniprot ID	Q15116
Description	PE-conjugated Anti-PD-1 (pembrolizumab biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

