

PRODUCT INFORMATION

Common Name	MPDL3280A,RG7446
Conjugate	PE-conjugated
Synonyms	PD-L1;CD274;B7-H1;PDCD1L1;PDCD1LG1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	PDL1
Uniprot ID	Q9NZQ7
Description	PE-conjugated Anti-PDL1 (atezolizumab biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

