

PRODUCT INFORMATION

Common Name	PF-06523435,hu24
Conjugate	PE-conjugated
Synonyms	CCK4
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 1 μ l/test
Formulation & Reconstitution	Liquid,PBS with,0.18%BSA□0.1%Proclin300
Host Species	Humanized
IgG type	Human IgG1(E356D,M358L) - kappa
Reactivity	Human
Target	PTK7
Uniprot ID	Q13308
Description	PE-conjugated Anti-PTK7(cofetuzumab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at 2°C-8°C for 6 months
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μ m) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 1µl/test PE-conjugated Anti-PTK7(cofetuzumab biosimilar) mAb on 293T cells (Blue histogram) or Raji cells (Red histogram).

