

PRODUCT INFORMATION

Common Name	PF-06523435,hu24
Conjugate	PE-conjugated
Synonyms	CCK4
Applications	Flow Cyt
Recommended Dilutions	1µl/test
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1(E239D, M241L)
Reactivity	Human
Target	PTK7
Uniprot ID	Q13308
Description	PE-conjugated Anti-PTK7(cofetuzumab biosimilar) mAb
Delivery	In Stock
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.






Figure 1. Flow cytometry analysis with 1µl/test PE-conjugated Anti-PTK7(cofetuzumab biosimilar) mAb on 293T cells (Blue histogram) or Raji cells (Red histogram).

