

PRODUCT INFORMATION

|                              |   |
|------------------------------|---|
| Common Name                  | BA302   |
| Conjugate                    | PE-conjugated   |
| Synonyms                     | NTRKR2  |
| Applications                 | Flow Cyt  |
| Recommended Dilutions        | Flow Cyt 1:100  |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA  |
| Host Species                 | Chimeric Humanized  |
| IgG type                     | IgG1(E239D, M241L)  |
| Reactivity                   | Human   |
| Target                       | ROR2  |
| Uniprot ID                   | Q01974  |
| Description                  | PE-conjugated Anti-ROR2(ozuriftamab biosimilar) mAb   |
| Delivery                     | Under Development   |
| Storage & Shipping           | Store at 2°C-8°C for 6 months   |
| Background                   | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage                        | Research use only   |

