

PRODUCT INFORMATION

| | |
|---|--|
| Common Name | SEG101, SelG1, crizanlizumab-tmca |
| Conjugate | PE-conjugated |
| Synonyms | P-Selectin;CD62P;SELP;GMP-140 |
| Applications | Flow Cyt |
| Endotoxin | Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization. |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Humanized |
| IgG type | Human IgG2 - Kappa |
| Reactivity | Human |
| Target | SELP |
| Uniprot ID | P16109 |
| Description | PE-conjugated Anti-SELP(crizanlizumab biosimilar) mAb |
| Delivery | Under Development |
| Storage&Shipping | Store at 2°C-8°C for 6 months |
| Sterility | Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μ m) prior to use. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

