Cat. No. DMC100445P



## **PRODUCT INFORMATION**

Clone ID **DMC445 Target** SLAMF1

**Synonyms** SLAMF1; CD150; CDw150; SLAM; IPO-3

**Host Species** 

PE-conjugated Anti-SLAMF1 antibody(DMC445); Description

IgG1 Chimeric mAb

**Delivery Under Development** 

**Uniprot ID** Q13291

Rabbit/Human Fc chimeric IgG1 IgG type

Clonality Monoclonal Reactivity Human **Applications** Flow Cyt

Recommended

Flow Cyt 1:100 **Dilutions** 

Purified from cell culture supernatant by affinity **Purification** 

chromatography

Formulation & Reconstitution

**DIMA Disclaimer** 

Liquid PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping Store at 2°C-8°C for 6 months

> SLAMF1 (Signaling Lymphocytic Activation Molecule Family Member 1) is a Protein Coding gene. Diseases associated with SLAMF1 include Measles and Subacute Sclerosing Panencephalitis. Among its related pathways are

**Background** Dendritic Cells Developmental Lineage Pathway

and T Cell Co-Signaling Pathway: Ligand-Receptor Interactions. Gene Ontology (GO) annotations related to this gene include antigen binding. An

important paralog of this gene is LY9.

**Usage** Research use only

PE-conjugated Conjugate

> All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are

> > Email: info@dimabio.com Website: www.dimabio.com

actively scrutinizing all patent application to

ensure no IP infringement.

