

PRODUCT INFORMATION

Common Name	XmAb-18087
Conjugate	PE-conjugated
Synonyms	SS-2-R;SS2-R;SS2R;SST2;SRIF-1
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 2μl/test
Formulation & Reconstitution	Liquid,PBS with,0.18%BSA□0.1%Proclin300
Host Species	Chimeric Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	SSTR2
Uniprot ID	P30874
Description	PE-conjugated Anti-SSTR2(tidutamab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at 2°C-8°C for 6 months
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μm) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 0.25µl/test PE-conjugated Anti-SSTR2(tidutamab biosimilar) mAb on Expi293 cells transfected with Human SSTR2 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

