

PRODUCT INFORMATION

| | |
|------------------------------|---|
| Common Name | MORAb-004 |
| Conjugate | PE-conjugated |
| Synonyms | CD248; CD164L1 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Chimeric Humanized |
| IgG type | IgG1(E239D,M241L,S286F) |
| Reactivity | Human |
| Target | TEM1 |
| Uniprot ID | Q9HCU0 |
| Description | PE-conjugated Anti-TEM1(ontuxizumab biosimilar) mAb |
| Delivery | Under Development |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

