

PRODUCT INFORMATION

Common Name	TSR-022
Conjugate	PE-conjugated
Synonyms	HAVCR2;TIM3;TIMD3;FLJ14428;KIM3
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG4 - Kappa
Reactivity	Human
Target	TIM3
Uniprot ID	Q8TDQ0
Description	PE-conjugated Anti-TIM3 (cobolimab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

