

PRODUCT INFORMATION

Common Name	MBG-453
Conjugate	PE-conjugated
Synonyms	HAVCR2, TIMD3
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG4
Reactivity	Human
Target	TIM3
Uniprot ID	Q8TDQ0
Description	PE-conjugated Anti-TIM3(sabatolimab biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

