

PRODUCT INFORMATION

Common Name	OPN-305
Conjugate	PE-conjugated
Synonyms	CD282;TIL4
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric
IgG type	IgG4
Reactivity	Human
Target	TLR2
Uniprot ID	O60603
Description	PE-conjugated Anti-TLR2(tomaralimab biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

