

PRODUCT INFORMATION

Common Name	PF 6480605,PF-6480605,RVT-3101,1D1-1.31,afimkibart
Conjugate	PE-conjugated
Synonyms	TL1;VEGI
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG1(L234A,L235A,G237A) - Kappa
Reactivity	Human
Target	TNFSF15
Uniprot ID	O95150
Description	PE-conjugated Anti-TNFSF15(PF-06480605 biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μ m) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

