

PRODUCT INFORMATION

Common Name	IMMU-132, hRS7-SN-38, hRS7-SN-38-ADC, hRS7-[CL-SN-38], hrS7-SN-38
Conjugate	PE-conjugated
Synonyms	TACSTD2;GA733-1;M1S1;TROP2
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 0.5μl/test
Formulation & Reconstitution	Liquid,PBS with,0.18%BSA□0.1%Proclin300
Host Species	Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	Trop2
Uniprot ID	P09758
Description	PE-conjugated Anti-Trop2 (sacituzumab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at 2°C-8°C for 6 months
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μm) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 0.5ul/test PE-conjugatedAnti-Trop2 (sacituzumab biosimilar) on T-47D cell line (Blue histogram) or 293T (Red histogram).

