

PRODUCT INFORMATION

Common Name	IMMU-132, hRS7-SN-38, hRS7-SN-38-ADC, hRS7-[CL-SN-38], hrS7-SN-38
Conjugate	PE-conjugated
Synonyms	TACSTD2;GA733-1;M1S1;TROP2
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 0.5µl/test
Formulation & Reconstitution	Liquid,PBS with,0.18%BSA□0.1%Proclin300
Host Species	Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	Trop2
Uniprot ID	P09758
Description	PE-conjugated Anti-Trop2 (sacituzumab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 0.5ul/test PE-conjugatedAnti-Trop2 (sacituzumab biosimilar) on T-47D cell line (Blue histogram) or 293T (Red histogram).

