

GPCR Flow Cytometry Antibodies Promotion List

Target	Cat.No.	Product Name	Applications
ADGRE2	DMC100370	Anti-ADGRE2 antibody(DMC370); IgG1 Chimeric mAb	Flow Cyt
C5AR2	DMC100482	Anti-GPR77 antibody(DMC482); IgG1 Chimeric mAb	Flow Cyt
CB1	DME100144	Anti-CB1 antibody(DM144); Rabbit mAb	ELISA; Flow Cyt
CCR1	DMC100465	Anti-CCR1 antibody(DMC465); IgG1 Chimeric mAb	Flow Cyt; WB
CCR2	DMC100447	Anti-CCR2 antibody(DMC447); IgG1 Chimeric mAb	Flow Cyt
CCR4	DMC100297	Anti-CCR4 antibody(DMC297); IgG1 Chimeric mAb	Flow Cyt
CCR5	DMC101003	Anti-CCR5 antibody(4F6); IgG1 Chimeric mAb	Flow Cyt
CCR6	DMC100477	Anti-CCR6 antibody(DMC477); IgG1 Chimeric mAb	Flow Cyt
CCR8	DME100180	Anti-CCR8 antibody(DM180); Rabbit mAb	ELISA; Flow Cyt
	DME100196	Anti-CCR8 antibody(DM196); Rabbit mAb	ELISA; Flow Cyt
CXCR1	DMC100470	Anti-CXCR1 antibody(DMC470); IgG1 Chimeric mAb	Flow Cyt
CXCR2	DMC100681	Anti-CXCR2 antibody(DMC681); IgG1 Chimeric mAb	Flow Cyt
CXCR3	DME100208	Anti-CXCR3 antibody(DM208); Rabbit mAb	ELISA; Flow Cyt
CXCR4	DMC100680	Anti-CXCR4 antibody(DMC680); IgG1 Chimeric mAb	Flow Cyt
CXCR5	DMC100679	Anti-CXCR5 antibody(DMC679); IgG1 Chimeric mAb	Flow Cyt; IHC
	DMC100679P	PE-conjugated Anti-CXCR5 antibody(DMC679); IgG1 Chimeric mAb	Flow Cyt
CXCR7	DMC100443	Anti-CXCR7 antibody(DMC443); IgG1 Chimeric mAb	Flow Cyt
FZD10	DMC100487	Anti-FZD10 antibody(DMC487); IgG1 Chimeric mAb	Flow Cyt
FZD4	DMC101111	Anti-FZD4 antibody(10E10); IgG1 Chimeric mAb	Flow Cyt
GIPR	DMC101001	Anti-GIPR antibody(3G5); IgG1 Chimeric mAb	Flow Cyt
GLP1R	DME101291	Anti-GLP1R antibody(1A4); Rabbit mAb	Flow Cyt
GPR20	DME101292	Anti-GPR20 antibody(1E10); Rabbit mAb	Flow Cyt
GPR75	DMC100368	Anti-GPR75 antibody(DMC368); IgG1 Chimeric mAb	Flow Cyt
	DMC100492	Anti-GPR75 antibody(DMC492); IgG1 Chimeric mAb	Flow Cyt
	DMC100368B	Biotinylated Anti-GPR75 antibody(DMC368); IgG1 Chimeric mAb	Flow Cyt
	DMC100492B	Biotinylated Anti-GPR75 antibody(DMC492); IgG1 Chimeric mAb	Flow Cyt
	DME101136	Anti-GPR75 antibody(12B6); Rabbit mAb	Flow Cyt
	DME101136B	Biotinylated Anti-GPR75 antibody(12B6); Rabbit mAb	Flow Cyt
	DMC100368P	PE-conjugated Anti-GPR75 antibody(DMC368); IgG1 Chimeric mAb	Flow Cyt
GPR87	DMC100478	Anti-GPR87 antibody(DMC478); IgG1 Chimeric mAb	Flow Cyt
	DMC100478B	Biotinylated Anti-GPR87 antibody(DMC478); IgG1 Chimeric mAb	Flow Cyt
GPRC5D	DME100091	Anti-GPRC5D antibody(DM91); Rabbit mAb	ELISA; Flow Cyt
	DME100091B	Biotinylated Anti-GPRC5D antibody(DM91); Rabbit mAb	ELISA; Flow Cyt
	DME100091P	PE-conjugated Anti-GPRC5D antibody(DM91); Rabbit mAb	Flow Cyt
GRPR	DME101290	Anti-GRPR antibody(3B12); Rabbit mAb	Flow Cyt
LGR4	DME101188	Anti-LGR4 antibody(1D9); Rabbit mAb	Flow Cyt
MC4R	DME101294	Anti-MC4R antibody(2H10); Rabbit mAb	Flow Cyt
SSTR2	DMC100428	Anti-SSTR2 antibody(DMC428); IgG1 Chimeric mAb	Flow Cyt
TSHR	DMC101228	Anti-TSHR antibody(1C8); IgG1 Chimeric mAb	Flow Cyt

 Email: info@dimabio.com

 Website: www.dimabio.com/www.dimabio.cn